

## Material required

- pH indicator
- BecKer
- Some tubes
- Water
- Vinegar
- HCI
- NaOH
- Ammonia
- Lemon
- Water with bicar





We make the substances and the indicator.

Vinegar HCl Lemon




Indicator becomes <u>red</u>





pH < 7



These are <u>acid</u> substances

Ammonia NaOH Water with Bicarbonate




Indicator becomes <u>blue</u>





pH > 7



These are <u>basic</u> substances Water




Indicator doesn't change color





pH = 7



This is a neutral substance



The indicator changes colour in each tests!